



CLIP RODDERS



RODDER

The Adept Lamb Rodder is used to separate the weasand (esophagus) away from the trachea and from other connecting tissue. The Adept Lamb Weasand Clip is then applied and pushed up into position using the Adept Lamb Rodder.

- > Designed for optimum separation of connecting tissue - no cutting required
- > Fast and easy to use
- > Lightweight and ergonomic
- > Polished finish for ease of cleaning and sterilization

Product Code	0142-050
Supplied	Single unit
Colour	White, Orange, Green (additional colours available on request)
Material	Stainless Steel / Plastic
Dimensions	278mm
Clips	Lamb Clip (0285-001,002,003,004) Lamb Clip Cart (0664-001,002,003,004)



WITH BLADE

- > Added Blade feature to allow quick trimming
- > Designed for optimum separation of connecting tissue - no cutting required
- > Lightweight and ergonomic
- > Polished finish for ease of cleaning and sterilization

Product Code	0142-090
Supplied	Single unit
Colour	White, Orange, Green (additional colours available on request)
Material	Stainless Steel / Plastic
Dimensions	278mm
Clips	Lamb Clip (0285-001,002,003,004) Lamb Clip Cart (0664-001,002,003,004)